
# TiMESH®

## pfm medical

pfm medical titanium gmbh

Postfach 2833 90013 Nuernberg Germany

Tel: +49 (0)911 215370-0

Fax: +49 (0)911 215370-30

www.pfmmedical.com

Federal law (USA) restricts this device to sale by or on the order of physician.

Patented by EP 0 897 997 A1, US 6,057,031, further patents pending

MK-0000726\_02

2011-02

## Instructions for use

### Description

Large-pore, non-absorbable prosthetic mesh made from monofilament polypropylene thread with surrounding titanized coating. The mesh is perfectly coordinated with the dynamics of body tissue in terms of resistance to tearing and elasticity. The mesh implants are available in various basis weights with defined shapes and dimensions.

### Indications for use

The TiMESH®, also known as TiMESH-TC®, along with Ethicon Inc. Prolene Soft Mesh, U.S.S.C.'s NonAbsorbable Polypropylene Surgical Mesh and Coated SURGIPRO Polypropylene Surgical Mesh are intended to be used for the reinforcement of tissue during surgical repair. Thus, the TiMESH® also known as TiMESH-TC® mesh and all of the predicates have the same intended use.

The TiMESH® is specifically indicated for laparoscopic and open surgery for the repair of direct or indirect inguinal, femoral, umbilical, and incisional hernias; rectal, vaginal, and apical prolapses; and parietal reinforcement of tissues and abdominal wall repair. The TiMESH® is specifically indicated for open surgery for the repair of direct or indirect inguinal and femoral hernias, parietal reinforcement of tissues, and abdominal wall repair. The TiMESH® has the same indications as a combination of the predicate devices. TiMESH® is a prescriptive device and should only be used by a licensed physician.

### Mode of action

The purpose of the implants is to strengthen tissue or bridge defects in soft tissue. Adsorption of a thin fibrous tissue layer, whose cells grow through the pores of the mesh fabric, results in the complete enclosure of the mesh material by the patient's own body tissue. Experimental investigations have shown that the implant remains soft and flexible in this context.

### Warnings

Complications which can occur during a reconstruction of the tissue walls with a wall reinforcement as a result of the foreign body reaction to plastics are greatly reduced by the titaniferous coating of TiMESH® implants. Although allergic reactions are significantly reduced, they cannot, however, be ruled out in principle (in particular with possible individual adaptation through cutting). We therefore recommend not damaging the titaniferous coating so as not to expose the plastic.

To guarantee the properties of the product, it must be ensured that sharp objects do not affect the mesh beyond the usual clinical application during the operation and that the titaniferous coating and plastic as carrier are not damaged.

**This product contains chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.**

### Safety precautions

- With incisional and abdominal wall hernias, the peritonealisation of the TiMESH® mesh implant at the end of the intervention should be as complete as possible.
- It is recommended to allow the hernia reinforcement to greatly overlap the edges of the hernial orifice.
- With a wide dissection and large tissue detachment, it is recommended to carefully check the haemostasis and avoid remaining cavities as far as possible, in order to facilitate the rapid integration of the hernia reinforcement.
- To counteract the inflammatory effect of foreign bodies, such as talcum, the implant should only be handled with unpowdered gloves.
- TiMESH® implants are supplied in sterile packaging. Carefully check whether the packaging is intact before use. On no account use the implant if the packaging is damaged, open or has been penetrated by moisture.
- The implant is sterile packed in a double package. The inner package should not be opened until the implant is used and then only handled with sterile gloves and sterile instruments.
- The TiMESH® implants are only sold to specialist surgeons. Implants are carried out entirely on the responsibility of the surgeon.
- Before using TiMESH® extralight hernia meshes, the physician should check the size of the hernial orifice and ensure a sufficient overlap.

### Contraindications

All the known usual contraindications and rules for soft tissue implants are valid for the use of TiMESH® tissue wall reinforcements.

These include, among others:

- Patients in the growth phase (the recognised ability of this reinforcement to extend can prove to be insufficient for the growth of children).
- Intervention in contaminated or inflamed locations.
- With large hernias the suitability of TiMESH® extralight mesh implants must be checked and a corresponding overlap of the hernial orifice ensured.
- The TiMESH® mesh implant is only intended for a single use and is not suitable for reuse. The properties of the product can change if resterilised.

### Application for the inguinal hernia indication

TiMESH® implants meet the demands of modern surgery and the current state of surgical mesh operation techniques, the foreign body reaction being clearly reduced as a result of the titaniferous layer in particular. They are used in both open (e.g. acc. to Lichtenstein etc.) and laparoscopic implantation (e.g. TAPP, TEP). Sufficient preparation, an adequate size of mesh and, if necessary, an anatomically suitable cutting to size must be ensured. To avoid a dislocation or distortion and to prevent the edges from rolling up, the mesh can be fixed with nonabsorbable threads or suitable clips in accordance with the respective operation technique, after placement free of tension, creases and cavities.

The excellent handling qualities of TiMESH® extralight hernia meshes ensure that rolling up or folding is not necessary before introduction into the trocar. The mesh is gripped in the middle with the forceps and introduced intracorporeally with a reduction sleeve.

With the extralight version trocars with a small diameter (≤ 5 mm) can be used. The mesh can be cut to size with scissors or a scalpel for individual requirements. Please observe our information on possible complications.

### Application for the incisional hernia indication

TiMESH® implants can be implanted either in the perineum between the abdominal wall and the peritoneum, or before the deep muscle layer. Can be used as a whole or individually adapted (please observe our information on possible complications).

The fasteners should be sufficiently close to each other, so that no physical weak point occurs during the surgical repair or at the beginning of a possible recurrence. The fixation should be at least 1 cm (3/8") from the edge of the mesh. Clips or suture can be used for fastening to the surrounding tissue.

### Storage

At room temperature in the original packaging away from light. Protect from moisture.

LOT	Batch number	🕒	Use by
STERILE	Sterilisation with ethylenoxide	⚠️	Follow instructions
REF	Order number	🚫	Do not reuse

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